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Stain and use	Hucker's		Carbol fuchsin		Kopeloff's	
	Reagent	Time	Reagent	Time	Reagent	Time
Initial stain	Crystal violet	30 s	Crystal violet	30 s	Alkaline crystal violet: flood with solution A; add 5 drops of solution B	2-3 min
Iodine	Gram's iodine	30 s	Gram's iodine	30 s	Kopeloff's iodine	≥2 min
Decolorizer	Acetone-alcohol	~1-5 s	95% ethanol	~30 s	3:7 acetone-alcohol: rinse immediately after applying	
Counterstain	Safranin*	30 s	Carbol fuchsin or 0.8% basic fuchsin	≥1 min	Kopeloff's safranin	10-30 s
Recommended use	General bacteriology		<i>Bacteroides</i> spp. <i>Fusobacterium</i> spp. <i>Legionella</i> spp. <i>Campylobacter</i> spp. <i>Brucella</i> spp. and other faintly staining gram-negative organisms		Anaerobes Diagnosis of bacterial vaginosis (Appendix 3.2.1-3)	

* Or, preferably, use 0.1 to 0.2% basic fuchsin as a counterstain (7).

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
Organism(s)	Gram stain morphology	Frequent sources	Comments and additional tests or media that may be included
<i>Actinomyces</i> spp.	Thin, beaded, branched gram-positive filaments; may be within sulfur granules with peripheral clubs	Cervicofacial, thoracic, abdominal, and pelvic abscesses and drainages; pleural fluid; bronchial washings	Modified acid-fast stain; if sulfur granules present, wash and crush in THIO
<i>Nocardia</i> spp.	Long, thin, branching, beaded, gram-positive or irregularly staining bacilli; in culture smears, may be pleomorphic with branching and coccoid forms	Sputum, bronchial washings, biopsy material, purulent exudates, CSF, blood	Modified acid-fast stain; if mixed microbiota, mycobacterial decontamination procedures may be used; plate incubated at 45°C may enhance recovery; use CYE ^a or Thayer-Martin agar.
<i>Mycobacterium</i> spp.	Gram-positive beaded or gram-neutral bacilli; often found inside macrophages; bacilli may be short to long, banded or beaded, and/or slightly curved; some species appear pleomorphic and coccoid	Respiratory tract, urine, blood, biopsy material, CSF	Confirm with acid-fast stain. Add AFB culture.
<i>Corynebacterium</i> spp.	Gram-positive pleomorphic, club shaped, irregularly staining bacilli or coccobacilli with palisading and/or angular arrangements	Blood, tissue aspirates, skin lesions, wounds, indwelling catheters, prosthetic heart valves, upper and lower respiratory tracts	Add selective and differential media for <i>Corynebacterium diphtheriae</i> , if suspected.
<i>C. jeikeium</i>	Often small coccobacilli resembling streptococci		
<i>Propionibacterium</i> spp.	Gram-positive, very pleomorphic "diphtheroid" forms that may branch	Blood, CSF, other body fluids, skin lesions	Common skin contaminant during needle aspiration
<i>Listeria monocytogenes</i>	Gram-positive small to medium coccobacilli that may be pleomorphic; occur in short chains or palisades; may be confused with corynebacteria or enterococci	CSF, blood, amniotic fluid, placental or fetal tissue	Direct wet mount for tumbling motility; if mixed microbiota, cold enrichment may be used
<i>Erysipelothrix rhusiopathiae</i>	In tissue, long, slender, gram-positive bacilli; in blood, small "coryneforms"	Skin lesions, biopsy material, tissue aspirates, blood	Associated with occupational or animal contact
<i>Lactobacillus</i> spp.	Medium, straight, uniform gram-positive bacilli with rounded ends; may form chains or spirals; sometimes short and coccobacillary	Usually involved in mixed infections; rarely from blood, CSF	Recovery may be improved by anaerobic incubation; normal vaginal, mouth, and gastrointestinal tract microbiota
<i>Bacillus</i> spp.	Medium to large square-ended gram-positive bacilli with parallel sides with or without spores; some species have spores that swell sides; may stain gram variable or gram negative with age	May be clinically relevant from any source in compromised patient or intravenous-drug abuser; also intraocular	Frequent culture contaminants; may cause ocular infections

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Organism(s)	Gram stain morphology	Frequent sources	Comments and additional tests or media that may be included
<i>Clostridium perfringens</i>	Gram-positive large "boxcars" with no spores; may stain gram negative	Blood, wounds, intra-abdominal	Add egg yolk agar incubated anaerobically; absence of inflammatory cells may indicate clostridial myonecrosis; normal gastrointestinal tract microbiota
<i>Clostridium</i> spp.	Gram-positive, -variable, or -negative bacilli with or without spores; bacilli may be large, slender and short, or long; sometimes form coils; often smaller than <i>Bacillus</i> spp.	Blood, intra-abdominal, wounds, abscesses	Normal gastrointestinal and genital tract microbiota
<i>S. pneumoniae</i>	Gram-positive cocci in pairs, lancet shapes, short chains	Lower respiratory tract, blood, CSF, sterile fluids	Quellung test may be used on selected clinical specimens.
<i>Enterococcus</i> spp.	Gram-positive cocci in pairs, short chains; may resemble pneumococci	Urine, wounds, blood, intra-abdominal abscesses	Normal gastrointestinal tract microbiota; common cause of superinfections in patients treated with expanded-spectrum cephalosporins
<i>Streptococcus</i> spp.	Round to oval gram-positive cocci, occasionally elongated; in pairs and/or short to long chains; nutritionally variant streptococci often seen as highly pleomorphic, gram-variable to gram-negative coccobacilli with pointed ends and spindle shapes	Blood, CSF, respiratory tract, multiple other sources	May be difficult to distinguish from corynebacteria and lactobacilli
<i>Aerococcus viridans</i>	Gram-positive cocci in pairs, tetrads, clusters	Blood, CSF	
<i>Staphylococcus</i> spp.	Gram-positive cocci in pairs, tetrads, clusters	Abscesses, drainages, wounds, respiratory tract, blood, tissue, sterile fluids, indwelling catheters	Normal microbiota, especially skin, nares
<i>S. aureus</i>	May often be characterized by very uniform, geometric clusters of small cocci, whereas coagulase-negative species are often irregular and more pleomorphic, with greater size variation		
<i>Rothia mucilaginosa</i>	Large gram-positive cocci in pairs, tetrads	Blood in compromised patients, peritoneal dialysates	Normal oral microbiota

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Organism(s)	Gram stain morphology	Frequent sources	Comments
<i>Enterobacteriaceae</i>	Straight thick bacilli; short to medium length with rounded ends; antimicrobial agent-affected organisms may be pleomorphic, filamentous, and/or irregularly staining	Urinary tract, multiple other sources	Includes organisms that cause gastroenteritis and bacterial dysentery; also normal gastrointestinal tract microbiota; nosocomial strains may be multiply resistant to antimicrobial agents
<i>Pseudomonas</i> spp.	Thin bacilli; medium length to long with rounded to pointy ends; often in pairs; antimicrobial agent-affected organisms may appear filamentous, coiled, and/or pleomorphic	Lower respiratory tract, wounds, eyes, multiple sites in compromised patients	Nosocomial strains may be multiply resistant to antimicrobial agents
<i>Haemophilus</i> spp., <i>Pasteurella</i> spp., fastidious gram-negative bacilli	Small coccoid to bacillary forms; pleomorphic; often with filamentous forms; may be faintly staining	Blood, sterile fluid (including CSF), lower respiratory tract, abscesses, wounds, eyes, genital tract	Inoculate CHOC.
<i>Legionella</i> spp.	Pleomorphic slender bacilli of variable lengths that may stain pale; may not take stain in clinical specimens	Lower respiratory tract	Add special growth media; direct fluorescent antibody stains and molecular probes available; acid wash method may be used to enhance recovery from specimens with mixed microbiota
<i>Fusobacterium nucleatum</i> , <i>Capnocytophaga</i> spp.	Long slender bacilli with tapered to pointed ends; "needlelike"; may be in pairs, end to end, or filamentous	Respiratory tract, wounds, blood, abscesses	Endogenous microbiota
<i>Fusobacterium necrophorum</i> , <i>Fusobacterium mortiferum</i> , or <i>Fusobacterium varium</i>	Highly pleomorphic bacilli with rounded to tapered ends; pale and irregularly staining, with bizarre forms and round bodies	Wounds, blood, abscesses	Endogenous microbiota
<i>Bacteroides</i> spp.	Pleomorphic straight bacilli with possible irregular to bipolar staining	Wounds, blood, abscesses	Direct fluorescent antibody stain available; endogenous microbiota
<i>Vibrio</i> spp.	Slightly curved to straight bacilli	Stool, wounds	If mixed microbiota, selective medium (TCBS) recommended
<i>Campylobacter</i> spp. (<i>Helicobacter</i> spp.)	Thin, curved bacilli, including S shapes, gull wings, and long spiral forms	Stool, blood, gastric biopsy samples	Microaerophilic or capnophilic atmosphere required; if mixed microbiota, 42°C incubation recommended for recovery of thermophilic species
<i>Acinetobacter</i> spp.	Medium to large cocci in pairs; occasionally coccoid, bacillary, and filamentous forms; often resistant to decolorization	Urine, lower respiratory tract, blood, sterile fluids, wounds, abscesses, tissues, stool	Multiple sites in compromised patients
<i>Neisseria</i> spp., <i>Moraxella catarrhalis</i>	Medium to large cocci in pairs and tetrads, coffee bean shaped; no bacilli seen	Genital tract, urine, lower respiratory tract, blood, sterile fluids, wounds, abscesses	If mixed microbiota, selective medium may be used to enhance recovery of <i>Neisseria gonorrhoeae</i>
<i>Veillonella</i> spp.	Tiny cocci in sheets or clumps	Wounds, blood	Endogenous microbiota

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Organism(s)	Gram stain morphology	Frequent sources	Comments
<i>Pneumocystis carinii</i>	Gram-negative spherical cysts (5–7 μ m) often containing rosette of eight gram-negative intracytic bodies, or cluster of gram-negative cysts surrounded by halos in background of amorphous gram-negative material	Open lung and transtracheal biopsy materials, bronchial washings and lavage specimens, sputum	Confirm with Grocott methenamine-silver, toluidine blue O, Gram-Weigert, or direct or indirect fluorescent antibody stain.
<i>Blastomyces dermatitidis</i>	Gram-variable, broad-based, thick-walled yeast cell with figure-eight appearance	Bronchial washings, sputum, purulent exudates, skin lesions	
Microsporidia	Gram-variable spherical cells (1–4 μ m), thicker at one end	Stool, respiratory, cornea, urine	Confirm with chromotrope stain
<i>Cryptococcus neoformans</i>	Partially or completely gram-positive round yeast cell with clear or red-orange halo; yeast cells may appear stippled or gram neutral	CSF, blood, biopsy material, sputum, bronchial washings, cutaneous lesions	Confirm capsule with India ink; direct antigen detection procedures may be used on CSF; inoculate urea slant.
<i>Candida</i> spp.	Gram-positive budding yeast cell with or without pseudohyphae; may also appear stippled or gram neutral	Sputum, urine, blood, biopsy material, vaginal discharge, upper respiratory tract	Endogenous microbiota
<i>Malassezia furfur</i>	Bottle-shaped yeast cells in compact clusters, usually with short hyphal elements	Skin scrapings, blood drawn through catheter lines, hyperalimentation fluids	Inoculate lipid-enriched medium.

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